

RADIOLOGICAL CHARACTERIZATION Laboratory Procedures

LORNA JEAN H. PALAD
Health Physics Research Unit
Philippine Nuclear Research Institute
Commonwealth Avenue, Quezon city
Philippines
3-7 December 2007

PURPOSE

- To measure activity concentration of materials taken from a shutdown nuclear research reactor
- To comply with regulatory requirements for clearance of materials for unrestricted or restricted reuse

RADIONUCLIDE INVENTORY

1. Neutron activated materials

- Located in the reactor structure and have been irradiated by neutrons
- Major activation products – ^3H , ^{14}C , ^{22}Na , ^{36}Cl , ^{39}Ar , ^{41}Ca , ^{54}Mn , ^{55}Fe , ^{59}Ni , ^{63}Ni , ^{60}Co , ^{65}Zn , ^{93}Mo , ^{93}Zr , ^{94}Nb , $^{108\text{m}}\text{Ag}$, $^{110\text{m}}\text{Ag}$, ^{125}Sb , ^{133}Ba , ^{134}Cs , ^{152}Eu , ^{154}Eu , ^{155}Eu , $^{166\text{m}}\text{Ho}$


2. Contaminated Materials

- typically found as surface contamination in reactor systems where fuel failure has occurred, on fuel handling systems, and in primary heat transport systems of any reactor
- Major fission products – ^{90}Sr , ^{99}Tc , ^{106}Ru , ^{129}I , ^{137}Cs , ^{144}Ce
- Major actinides - ^{238}Pu , ^{239}Pu , ^{241}Pu , ^{241}Am , ^{242}Cm , ^{244}Cm , ^{232}U , ^{233}U , ^{234}U , ^{235}U , ^{236}U , ^{238}U

MATERIALS TO BE SAMPLED

- Material samples
 - concrete, steel, resins, wood, tiles, paint chips, foam
 - pipes, sludge, pool water, etc...
- Environmental samples
 - soil/sediment, grass, water
- Biological
 - urine, excrete


Fig 1. Analytical Process



Field Sampling

Sampling design should be carefully planned

Sampling is the source of biggest contribution to error in analysis



Sample Receipt and Inspection

Sample Receipt

- receive in designated areas
- Check accompanying documents
- Conduct radiological survey

Sample Inspection

- Check physical integrity of package and container
- Confirm sample identity
- Confirm field preservation, if any

Laboratory Sample Tracking

- Sample log-in
- Sample tracking during analysis
- Storage of sample

Detection Methods, MDAs and Release Criteria for Radionuclides

Isotope	Half-life	Emission	Detection Methods	MDA, Bq/g	Clearance Level, Bq/g
H-3	1.2E+01	β^-	Liquid Scintillation	10	100
C-14	5.7E+03	β^-	Liquid Scintillation	1	1
Cl-36	3.0E+05	β^-	Liquid Scintillation	1	1
Ca-41	1.0E+05	EC	Liquid Scintillation	(1 to 10)	
Mn-54	8.6E-01	EC, γ	Gamma spectrometry	0.5	
Fe-55	2.7E+00	EC, X	Liquid Scintillation or X-ray spectrometry	10	
Co-60	5.3E+00	β^- , γ	Gamma spectrometry	0.5	0.1
Ni-59	7.5E+04	EC, X	X-ray spectrometry	10	100
Ce-144	8.0E-01	β^- , γ	Gamma spectrometry	0.5	
Ni-63	1.0E+02	β^-	Liquid Scintillation	1	100
Cs-134	2.1E+00	β^- , γ	Gamma spectrometry	0.5	0.1
Cs-137	3.0E+01	β^- , γ	Gamma spectrometry	0.5	0.1

Isotope	Half-life	Emission	Detection Methods	MDA, Bq/g	Clearance Level, Bq/g
I-129	1.6E+07	β^-	ICPMS (or X ray spectrometry)	0.007	0.01
Sr-90	2.9E+01	β^-	Beta counting or Liquid scintillation	1	1
Zr-93	1.5E+06	β^-	ICPMS	0.1	
Nb-93m	1.6E+01	Γ, X	X ray spectrometry or liquid scintillation	10	
Nb-94	2.0E+04	β^-, γ	Gamma spectrometry (or ICPMS)	0.5 (7)	0.1
Mo-93	3.5E+03	EC,X	X ray spectrometry or liquid scintillation	10	10
Tc-99	2.1E+05	β^-	ICPMS	0.6	
Ru-106	1.0E+00	β^-, γ	Gamma spectrometry	0.5	
Ag108m	1.3E+02	EC, γ	Gamma spectrometry	0.5	
Ag-110m	7.0E-01	β^-, γ	Gamma spectrometry	0.5	
Sb-125	2.8E+00	β^-, γ	Gamma spectrometry	0.5	
Ba-133	1.1E+01	EC, X, γ	Gamma spectrometry	0.5	

Isotope	Half-life	Emission	Detection Methods	MDA, Bq/g	Clearance Level, Bq/g
Eu-152	1.3E+01	EC, β^-, X, γ	Gamma spectrometry	0.5	0.1
Eu-154	8.6E+00	β^-, X, γ	Gamma spectrometry	0.5	0.1
Eu-155	4.8E+00	β^-, X, γ	Gamma spectrometry	0.5	1
Ho-166m	1.2E+03	β^-, X, γ	Gamma spectrometry	0.5	
U-234	2.5E+05	α, X	Alpha spectrometry	0.02	
U-235	7.0E+08	α, γ	ICPMS	0.0001	
U-238	4.5E+09	α	ICPMS	0.00001	
Pu-238	8.8E+01	α, X	Alpha spectrometry	0.02	
Pu-239	2.4E+04	α	Alpha spectrometry	0.02	0.1
Pu-241	1.4E+01	β^-	Liquid scintillation	1	
Am-241	4.3E+02	α, X, γ	Alpha spectrometry	0.02	
Cm-242	4.5E+01	α, X	Alpha spectrometry	0.02	
Cm-244	1.8E+01	α, X	Alpha spectrometry	0.02	

Gamma Spectrometry

- Nondestructive analysis of gamma emitting radionuclides
- Requires minimum sample preparation
- Measure gamma radiation using scintillation detectors (NaI (Tl)) crystals and semiconductor detectors (High purity germanium detectors)

I. Sample Preparation

1. Liquid Samples

- Evaporation
- Adsorption with AMP (Cs)/Precipitation with MnO₂ (other metals except Cr)
- Ion-exchange resin

2. Solid Samples

- Dry, grind, sieve, mix to homogenize, get sub sample

3. Air particulate filters, Swipe/Smear Samples

- No pretreatment
- Leach with dilute HNO₃ or 0.05M EDTA

- II. Standard Calibration
- III. Spectrometer Calibration
 1. Energy Calibration
 2. Efficiency Calibration
 3. Quality Control Check
- IV. Sample Measurement
 1. Count a blank sample
 2. Determine minimum detectable activity (MDA)
 3. Count sample
- V. Statistical Treatment of Data
 1. Calculate Activity including correction factors
 2. Calculate total uncertainty

Alpha Spectrometry

- Analysis of alpha emitting radionuclides
- Requires radiochemical separation
- Use of Si semiconductor detectors

I. Sample Treatment

1. Liquid Samples

- Acidify and filter
- May reduce volume by evaporation or ion exchange

2. Solid Samples

- Drying, grinding, sieving, homogenize
- May reduce volume by ashing
- Sample decomposed with concentrated acids, e.g. HCl, HNO₃, HF; H₂O₂ or alkali fusion with K₂CO₃
- Sample can be leached with dilute nitric acids or complexing agents

3. Air Particulate or Swipe/Smear Samples

- Cellulose and membrane filters can be dissolved in mineral acids
- Glass-fiber filters can be leached with dilute acids or chelating agents

I. Preconcentration

- Precipitation


II. Chemical Separation

- Anion exchange

III. Source Preparation

- Electrodeposition
- Precipitation with cerium or lanthanum

IV. Counting



Liquid Scintillation Counting

- Gross alpha/beta analysis of samples
- Analysis of beta emitters
- Analysis of alpha emitters
- May require radiochemical separation

I. Sample Preparation

1. Liquid Samples

- Distillation/electrolytic enrichment (^3H)
- Oxidation/ CO_2 trapping (^{14}C)
- Reduce volume by evaporation, ion exchange, precipitation

2. Solid Samples

- Dissolution with concentrated mineral acids (for activated and geogenic samples)
- Leaching with dilute mineral acids (for surface contaminated samples)
- Reduce volume by ashing at high temperatures (except volatile radionuclides, e.g. ^{129}I , ^{99}Tc)

3. Air Particulate/Swipe/Smear Samples

- Cellulose and membrane filters can be dissolved in concentrated mineral acids
- Glass fiber filters can be leached with dilute mineral acids or chelating agents
- Filter paper can be placed directly in LSC vial with addition of liquid scintillant

II. Setting of Measurement Conditions for Radionuclide

- Determination of Figure of Merit (FOM)

III. Drawing of Quench Correction Curve


- Preparation of quench standards
- Draw quench correction curve
- Efficiency determination

IV. Counting of Sample

V. Calculation of Activity

- Calculate Average, $X \pm 3\sigma$
- Correct activity using decay correction factor

Tritium Analysis by Liquid Scintillation Counting



THANK YOU